

TOC 1500

Gesamtorganischer Kohlenstoffgenerator

Artikelnummer: 42-1030




Beschreibung

Der TOC1500HP wurde für eine zuverlässige Luftversorgung von TOC-Analysatoren entwickelt und ist ein robustes System, das in Labors zur Erzeugung von sauberem, trockenem CO₂ und CH₄-frei eingesetzt wird



Generator Details

- **Gasarten:** Null-Luft
- **Integrierter Kompressor:** Nein



Hauptmerkmale

- Produziert saubere, trockene, CO₂- und CH₄-freie Luft
- Ideal für NPOC TOC-Analysatoren
- Funktioniert bei Bedarf rund um die Uhr. Wandmontage für kleine Labore
- Wenige bewegliche Teile bedeuten wenig Wartung und gewährleisten eine lange Lebensdauer des Generators
- Im Betrieb völlig geräuschlos
- Minimale Einrichtung erforderlich
- 12 Monate umfassende Vor-Ort-Garantie

Technische Spezifikationen

TOC 1500 Gesamtorganischer Kohlenstoffgenerator

Gasart	Null-Luft
Maximaler Gasfluss:	1500cc/min
Max. Ausgangsdruck:	80psi/5,5bar
Max. Ausgangsdruck:	80psi/5,5bar
Maximale Reinheit:	<1 ppm CO ₂ / 2.0ppm THC
Minimaler Eingangsdruck:	90psi/6,2bar
Maximaler Eingangsdruck:	145psi/10bar
Anforderungen an die Wasserreinheit:	N/A
Wasserverbrauch:	TBC
Stromspannung:	100-230 ± 10% VAC
Frequenz:	50 / 60 Hz
Strom:	8 Amps
Akkreditierungen:	CE
Größe (HxBxT) mm	660 x 250 x 160 mm
Größe (HxBxT) Zoll	26 x 9,8 x 6,3 inches
Gewicht des Generators	12kg / 26,4lbs

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

