

Corona Air Compressor

Funktioniert mit Corona 1010
Stickstoffgenerator

Artikelnummer: 85-0555 (230v) | 85-
1555 (110v)



Beschreibung

Der Druckluftkompressor Corona Air Compressor wurde so konstruiert, dass er sauber auf einen Corona 1010 Stickstoffgenerator aufgesetzt werden kann, falls das Labor keinen Zugang zu externer Luftzufuhr hat.




Generator Details

- **Gasarten:** Null-Luft
- **Integrierter Kompressor:** Ja

Hauptmerkmale

- Mit Corona 1010-Stickstoff stapelbar Bietet Luft für das Corona-System für Labore ohne externe Versorgung
- Gas wird auf Anfrage geliefert, sodass der Generator nach Ihrem Zeitplan arbeitet. 12 Monate umfassende Vor-Ort-Garantie



Technische Spezifikationen

Corona Air Compressor

Gasart	Null-Luft
Maximaler Gasfluss:	30L/min
Max. Ausgangsdruck:	110psi/7,6bar
Max. Ausgangsdruck:	110psi/7,6bar
Minimaler Eingangsdruck:	0
Maximaler Eingangsdruck:	110psi/7,6bar
Energieverbrauch:	572 VA / 736 VA
Stromspannung:	110 / 230 V
Frequenz:	50 / 60 Hz
Strom:	5.2 A / 3.2 A
Wärmeabgabe:	1950 VA / 2510 VA
Maximale Betriebstemperatur:	35°C / 95°F
Partikel:	<0.1µm
Akkreditierungen:	CE
Größe (HxBxT) mm	406 x 417 x 540 mm
Größe (HxBxT) Zoll	16 x 16,4 x 21,3 inches
Gewicht des Generators	42kg / 92,4lbs

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

