

# Precision Zero Air 7L


Flammenunterstützungsgenerator

Artikelnummer: 60-0008 (110V) | 60-0007 (230V)



## Beschreibung

Die Precision Zero Air 7L n- Nullluftgeneratoren sind speziell sauber, trocken, kohlenwasserstofffreier Luft zu versorgen als Flammstützgas für GC werden verwendet, um sowohl bei Standard- und Nachweisgrenzen verfolgen.




## Generator Details

- **Gasarten:** Null-Luft
- **Integrierter Kompressor:** Nein

## Hauptmerkmale

- Integriert in vorhandene Precision-Produkte mit demselben modularen, stapelbaren Design
- Kohlenwasserstoffgehalt (als Methan) & lt; 0,05 ppm für marktführende Reinheit
- Bis zu 7 l / min Zero Air
- Bis zu 5,5 bar Ausgangsdruck
- Mindestwartungsanforderungen, kein teurer Austausch der Katalysatorkammer
- Stark reduziertes Risiko, dass Verunreinigungen in das System gelangen, und kein Risiko, dass das Gas ausgeht
- Stilvolle und zeitgemäße Ästhetik und funktionales Design
- Statusanzeige-Beleuchtungsfunktion (bernsteinfarbenes Aufwärmen, grün bereit)
- Kompatibel mit Präzisionsluftkompressor
- 12 Monate umfassende Vor-Ort-Garantie



## Technische Spezifikationen

### Precision Zero Air 7L Flammenunterstützungsgenerator

Gasart	Null-Luft
Maximaler Gasfluss:	7L/min
Max. Ausgangsdruck:	80psi/5,5bar
Kohlenwasserstoffgehalt:	<0.05ppm
Max. Ausgangsdruck:	80psi/5,5bar
Gasanschluss:	1 x 1/4" BSPP
Minimaler Eingangsdruck:	110psi/7,6bar
Maximaler Eingangsdruck:	125psi/8,6bar
Min Air Inlet Flow:	7 L/min
Minimale Zuluftqualität:	ISO 8753-1:2010 Class 1.4.1 or Precision Air Compressor
Startzeit:	60 mins
Energieverbrauch:	594 W
Stromspannung:	110V / 230V
Frequenz:	50 / 60 Hz
Strom:	3.3 - 5.4 A
Wärmeabgabe:	2590
Maximale Betriebstemperatur:	35°C / 95°F
Akkreditierungen:	CE, FCC, CSA
Größe (HxBxT) mm	256 x 380 x 540 mm
Größe (HxBxT) Zoll	10,1 x 15 x 21,3 inches
Gewicht des Generators	25kg / 55lbs

## [PEAK Protected]™

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001

